


New LED Pear Lead

news



New Intercall LED Pear Leads give residents and users added peace of mind and reassurance

Now Fitted with re-assurance LED

The latest update to the market-leading Intercall Pear Lead features a built-in call re-assurance LED, which gives the user or resident an immediate visual indication that their call has successfully been generated. In addition, the button has been re-designed with a new crisp design and wipe clean finish.

The new series of lead builds on the success of earlier models – but gives the user added peace of mind and reassurance that help is on its way. Once the call has been made, the reassurance LED will not switch off until it has been reset at the call point.

Easy to fit and simple to use - The leads have been designed to clip onto the resident's clothing or bedding for immediate availability and ease of use, whilst the other end of the lead is plugged into the jack point in the room. In addition, the right-angled jack plug reduces the chances of accidental damage or removal by the movement of furniture. In the event that the lead is pulled out of its socket – a standard alarm is automatically generated.

Greater flexibility - For additional flexibility, the new Pear Leads are available in a choice of four different lengths: 2m, 4m, 6m and 8m, so residents have complete freedom to move around the room easily without the lead becoming detached.

Attractively styled and extremely simple to use, the new Intercall Pear Leads enable residents or patients with virtually all levels of disability to use the system with ease.

For more information on the new LED series or Pear Leads, call 0500 212123, or visit the Intercall web site at www.intercall.co.uk



Tel: 01403 713121

E-mail: sales@lismore.uk.com

www.intercall.co.uk

PearLead