

**NEW LIMPROG USB**

# product update



**Gives you even greater flexibility and control for your entire INTERCALL system**

## Limprog USB

The all new Intercall LIMPROG USB interface system is designed to give managers of nursing homes and hospitals even greater control over their Intercall systems by providing laptop control and management of the entire nurse call system.

The system, which replaces the earlier serial port unit, simply plugs into a USB port and enables remote programming and direct input for changing data such as patient names, care notes and other important information which can be shared by the entire system.

### CLEAR, EASY TO USE INTERFACE -

The Windows™ based software will run on any PC and enables you to process full system scans and download all system information such as number of devices connected, patient records and detailed site information including service and maintenance records.

In addition, the new Limprog USB integrates fully and easily with your operating system to give a sharp appearance and customisable look and feel.

The system also includes a network monitoring facility which can record all system activity together with event dates and times, which can be viewed in the form of an interactive 'virtual map' with extensive filtering which enables engineers to diagnose faults quickly and effectively.

Furthermore, as a special introductory offer, if you return your old style Lim, we will upgrade it for a special reduced price. To find out more about the new Limprog USB, simply call us FREE on 0500 212123, e-mail us at [sales@intercall.co.uk](mailto:sales@intercall.co.uk).



**Intercall Nursecall Systems**

Tel: 0500 212123

E-mail: [sales@lismore.uk.com](mailto:sales@lismore.uk.com)